



*Innovative Cell Technologies, Inc.*  
6790 Top Gun St. Suite 1, San Diego, CA 92121  
T.858 587-1716 F. 858 453-2117 E-mail.  
info@innovativecelltech

## Protocol for the use of ACCUMAX™ in Primary Tissue Dissociation

This protocol for using ACCUMAX™ to dissociate cells from primary tissue is a general-purpose protocol and may not be applicable to all tissue types. The individual/investigator needs to optimize the conditions for his/her tissue specimens. Keep in mind that ACCUMAX™ is a powerful enzyme mixture that can potentially dissolve not only the connective tissue of solid tissue but some fragile cell types as well if not closely monitored.

### MATERIALS

#### *Sterile:*

ACCUMAX™ (Should be defrosted overnight in the refrigerator or in a bucket of room temperature water - *not a 37°C bath*)  
DPBS (calcium and magnesium free)  
Culture medium, i.e., DMEM/F12 with 10 – 20% FBS (or other appropriate media)  
Pipettes - 1 ml, 10 ml  
Petri dishes -100 mm, non-tissue culture grade  
T25 culture flasks  
Centrifuge tubes, 15-50 ml, depending upon the amount of tissue being processed  
Scalpels  
Forceps

#### *Non-sterile:*

Platform rocker  
Trypan Blue  
Microscope  
Centrifuge

### PROCEDURE:

1. Transfer the tissue to a petri dish containing fresh, sterile DPBS, and rinse.
2. Transfer the tissue to a second dish; dissect off unwanted tissue, such as fat or necrotic material.
3. Using two crossed scalpels or a scalpel and forceps, cut the tissue into small pieces approximately 1 mm in size.
4. Transfer the tissue pieces to a 15 or 50 ml sterile centrifuge tube containing fresh, sterile DPBS.
5. Allow the pieces to settle and carefully remove the supernatant. Repeat this wash step two times.
6. Transfer the tissue pieces to a fresh petri dish and add enough ACCUMAX™ to the plate to cover tissue.
7. Incubate the samples on a platform rocker **at room temperature** 5 to 60 minutes. The tissue will “smear” on the bottom of the dish when the disaggregation is effective.
  - To release more cells, gently agitate the sample by pipetting several times.
  - It is best to check cell viability several times during the incubation using Trypan blue.
8. Once disaggregation is complete, transfer the cells to a sterile centrifuge tube and centrifuge at 300 x g to pellet the cells and to remove the cell debris if desired.
9. Carefully remove the supernatant and re-suspend the cell pellet in 5 ml of DMEM/F12 containing 10 – 20% FBS (or other appropriate media). Seed in a T25 flask. Replace the media after 48 hours.



*Innovative Cell Technologies, Inc.*  
6790 Top Gun St. Suite 1, San Diego, CA 92121  
T.858 587-1716 F. 858 453-2117 E-mail.  
info@innovativecelltech

## **Primary Tissue Dissociation Protocol, cont. page 2**

### **ALTERNATIVELY**

If cell isolation is from a soft tissue (such as liver):

1. Transfer the tissue to a petri dish containing fresh, sterile DPBS, and rinse.
2. Transfer the tissue to a second dish; dissect off unwanted tissue, such as fat or necrotic material. Add 1 – 2 ml of ACCUMAX™ and use forceps to gently “tease” the cells into the ACCUMAX™.
3. Residual connective tissue may be separated by allowing the pieces to settle or by filtration, if desired.
4. Centrifuge the sample at 300 x g to pellet the cells and to remove cell debris if desired.
5. Carefully remove the supernatant and re-suspend the cell pellet in 5 ml of DMEM/F12 containing 10 – 20% FBS (or other appropriate media). Seed in a T25 flask. Replace the media after 48 hours.